

Wolfi's law

General Use Protocol of the iMRS and the MRS2000+ designo®


Full Body Mat:

Use the full body mat two to three times a day for 8 minutes.

(Note: The cellular resonance lasts 4 to 6 hours)

- ✓ In the Morning set the level at 25
- ✓ In the Evening set the level at 10

Settings for the Probe or Pillow
(2-3 times daily)



Chronic
(older)

Acute
(new)

10 - 25
8 min

25 - 50
16 - 24 min

25 - 50
8 min

50 - 100
16 - 24 min

50 - 100
8 min

100 - 400
16 - 24 min

These setting are for:
very chronic conditions, pains
or symptoms that have been
bothering you for a long time.

These setting are for:
acute; new inflammatory
conditions, pains or symp-
toms. Situations that are
relatively new.

Very old persons
Very sick persons
Very fragile persons
etc...

Very inflamed areas
Very acute or severe pain
etc...

Speak to your representative if you have any questions. They are there to assist you with any questions. Do not adjust settings on machine without speaking to your representative first.

Technician Name:

Technician Number:

Important Note:

drink 6-8 glasses of water a day. This increases effectiveness of MRS sessions. (or half your body weight in ounces)

